DNA and RNA Extraction

General Information

We are pleased to offer a wide array of nucleic acid extraction options. Please note that the ones listed below reflect our most requested options but are not limited to these sample types. Please feel free to reach out to our project managers (pm.genome@mcgill.ca) in order to discuss your specific needs.

DNA extraction offered for:

- Whole blood 400 ml
- Whole blood 4-5 ml
- Buffy coat 200 µl
- Saliva 2 ml
- FTA cards
- FFPE
- Frozen tissue 10-30 mg/ sample
- Cell culture cells

RNA extraction offered for:

- Whole blood in Pax gene tubes
- Frozen cell pellets
- Frozen cell lysates
- Frozen tissue 10-30 mg/ sample
- Tissue or cells in RNAlater

Acceptable formats for each sample type

DNA extraction	Accepted container type
Whole blood 400 ml	2D Barcoded tubes (Matrix, FluidX, Eppendorf) Perkin Elmer 2ml Riplate (96) Deep well plates
Whole blood 4-5 ml	Vacutainer with anticoagulant; EDTA
Buffy coat 200 µl	2D Barcoded tubes (Matrix, FluidX, Eppendorf) Perkin Elmer 2ml Riplate (96) Deep well plates
Saliva 2 ml	Oragene•DNA OG-600 from DNA Genotek Oragene•DNA OG-500 from DNA Genotek
FTA cards	Whatman FTA Cards Punchouts (2-3x 5mm) from FTA card in 1.5 ml DNALoBind tubes
FFPE	2-3x 50µm microtome cut FFPE slices in 1.5 ml DNALoBind tubes
Frozen tissue	10-30mg in 1.5 ml DNALoBind tubes (PREFERRED) 1.5 - 2.0 ml screw cap tubes
Cell culture cells	Perkin Elmer 2ml Riplate (96) Deep well plates 2D Barcoded tubes (Matrix, FluidX, Eppendorf) 1.5 ml DNALoBind tubes 1.5 - 2.0 ml screw cap tubes

RNA extraction	Accepted container type
Whole blood in Pax gene tubes	PAXgene RNA Tubes for 2.5 ml whole blood
Frozen cell pellets	Perkin Elmer 2ml Riplate (96) Deep well plates 2D Barcoded tubes (Matrix, FluidX, Eppendorf) 1.5 ml DNALoBind tubes 1.5 - 2.0 ml screw cap tubes
Frozen cell lysates	Perkin Elmer 2ml Riplate (96) Deep well plates 2D Barcoded tubes (Matrix, FluidX, Eppendorf) 1.5 ml DNALoBind tubes 1.5 - 2.0 ml screw cap tubes
Frozen tissue 10-30 mg/sample	10-30 mg in 1.5 ml DNALoBind tubes (PREFERRED) 1.5 - 2.0 ml screw cap tubes
Tissue or cells in RNAlater	10-30 mg in 1.5 ml DNALoBind tubes (PREFERRED) 1.5 - 2.0 ml screw cap tubes

Useful links for accepted formats

1. 1.5 or 2.0 ml screw cap tubes, <https://www.fishersci.ca/ca/en/catalog/search/products?keyword=screw+cap+tubes>
2. Thermo Scientific™ Matrix™ 2D Barcoded Storage Tubes. Offered in three sizes, 0.5mL, 0.75mL, and 1.4mL, <http://www.thermofisher.com/order/catalog/product/3710>
3. FluidX 96-Well Format Sample Storage Tubes with Screw Cap and 2D Barcode External thread. Offered in three sizes, 0.3mL, 0.7mL, and 1.0mL, <http://www.fluidx.eu/96-well-format-sample-storage-tubes-with-2d-barcode.html>
4. Eppendorf 2D barcode Cryovials, <https://online-shop.eppendorf.ca/CA-en/Freezers-44537/Sample-Management-376575/Eppendorf-CryoStorage-Vials-PF-758909.html>
5. Perkin Elmer 2ml Riplate Deep 96 well plates Cat#43001-0120, http://www.analisisavanzados.com/modules/mod_cat/Ritter%20Robotic%20Puntas%20conductivas%20para%20Robots.pdf
6. Oragene collection tubes, <https://www.dnagenotek.com/ROW/industries/academic-research-human-genomics.html>
7. FTA cards, <https://www.sigmaaldrich.com/labware/labware-products.html?TablePage=108827546>
8. PAXgene RNA Tubes, <https://www.bdbiosciences.com/us/applications/blood-collection/cell-biomarker-preservation/paxgenereg-blood-rna-tube/p/762165>
9. Other formats will be discussed in a case-by-case basis.

Sample Submission Manifest

All samples must be accompanied by a sample submission sheet. Please follow the instructions within the file and fill in all appropriate fields. (Figure 1.). You can contact us if you need to provide more specifications than the sheet allows.

This completed document must be returned to us by email at the same addresses as for the tracking number and a hard copy must be included in your package.

Preparation for Shipping Samples

- Frozen samples must be shipped on dry ice or in a vapour phase cryoshipper in labeled boxes. Please ship no later than Tuesday morning in order to assure that samples arrive promptly and are still well frozen. Collaborators will be notified if samples arrived unthawed or show any sign of leakage or damage.
- Please follow EHS guidelines for shipping of dangerous goods including dry ice, biohazardous materials

McGill University Genome Center,
740 Docteur Penfield Avenue, Room 2304,
H3A 0G1, Montreal, QC.
c/o Corinne Darmond or Ariane Boisclair

Special Consideration for Level 2 Biological Samples

We are pleased to offer extractions on materials which require a Level 2 laboratory. Our BCL2+ lab comes equipped with the necessary equipment to handle samples which cannot be handled in our standard processing facility. We work closely with McGill’s Environmental Health and Safety team to ensure that every new potential pathogen is dealt with in a manner which is both safe for our personal as well as optimal for the samples themselves.

When shipping Level 2 biological samples, please refer to **the EHS TDG document** in compliance with provincial and federal regulations.

QC of Samples Post Extraction

DNA: Extracted DNA samples will be quantified using a PicoGreen assay if the number of samples extracted exceeds 24 samples. For 24 samples and under a Qubit quantification assay will be performed. Absorbance values (A260/A280 and A230/260) are not part of our regular QC but can be requested to the project managers.

Extracted DNA samples will be run on a 1% agarose gel to ensure that the extracted sample is not degraded. If a more accurate analysis of fragment size or DIN /RIN is requested, the samples can be run either on the Agilent TapeStation, Perkin Elmer LabchipGXT or the Agilent FemtoPulse. Please make sure to include this request during the discussion with our project managers.

RNA: For extracted RNA samples, QC will be performed using the Labchip GXT RNA kit for both quantification and assessment of quality (RIN) if the number of extracted samples exceeds 16. If the number of extracted samples is 16 or less, quantification will be performed using Qubit RNA assays followed by a TapeStation assay for evaluation of RNA quality (RIN).

Figure 1. Sample submission sheet

Sample Submission Template

Heading Rules
Only use the followin characters for Sample name and Barcode: a-z, A-Z, 0-9, period (.), dash (-), underscore (_), .

Biospecimen Type	Sample Name	Alias	Cohort	Experimental Group	Taxon	Container Kind	Container Name	Container Barcode	Location Barcode	Location Coord	Individual ID	Sex	Pedigree	Mother ID	Father ID	Volume (uL)	Conc. (ng/uL)	Collection Site	Tissue Source	Reception Date	Phenotype	Comment	

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Sample Submission Template

Naming Rules
Only use the follow characters for Sample name and Barcode: a-z, 0-9, period (.), dash (-), underscore (_)

Biospecimen Type	Sample Name	Alias	Cohort	Experimental Group	Taxon	Container Kind	Container Name	Container Barcode	Location Barcode	Location Coord	Individual ID	Sex	Pedigree	Mother ID	Father ID	Volume (uL)	Conc. (ng/uL)	Collection Site	Tissue Source	Reception Date	Phenotype	Comment